

Meetpunt	I	I Wide	II	II Wide	III	III Wide	IV	IV Wide	V	V Wide	VI	VI Wide
cT (AT Maxi)	> 80 cm		> 90 cm		> 90 cm		> 110 cm		> 110 cm		> 120 cm	
cH (AT Maxi)	> 95 cm		> 105 cm		> 105 cm		> 125 cm		> 125 cm		> 135 cm	
cG siliconenband	43-48	49-57	45-52	53-62	49-56	57-67	53-60	61-72	56-64	65-77	60-68	69-81
cG	43 - 57		45 - 62		49 - 67		53 - 72		56 - 77		60 - 81	
cF	39 - 52		41 - 56		44 - 60		47 - 65		50 - 69		53 - 73	
cE	30 - 37		33 - 40		35 - 43		37 - 45		39 - 48		41 - 51	
cD	27 - 33		29 - 36		32 - 39		34 - 42		36 - 45		38 - 48	
cC	28 - 34		30 - 37		33 - 40		35 - 43		37 - 46		39 - 49	
cB1	23 - 27		24 - 29		26 - 32		29 - 35		31 - 37		33 - 39	
cB	18 - 20		20 - 22		22 - 24		24 - 26		26 - 28		28 - 30	
cY	26 - 31		28 - 33		29 - 35		31 - 37		32 - 38		33 - 40	
cA	17 - 22		19 - 24		21 - 26		23 - 29		25 - 32		27 - 34	

Lengte D

ℓD petite: 34 - 38 cm
ℓD regular: 38 - 45 cm

Lengte G


ℓG petite: 60 - 69 cm
ℓG regular: 69 - 80 cm
ℓG long*: 80 - 90 cm

Lengte K (voor panty)

ℓK petite: 65 - 75 cm
ℓK regular: 75 - 85 cm

Maten zijn in cm

*Alleen beschikbaar in JOBST Classic



Point de mesure	I	I Wide	II	II Wide	III	III Wide	IV	IV Wide	V	V Wide	VI	VI Wide
cT (AT Maxi)	> 80 cm		> 90 cm		> 90 cm		> 110 cm		> 110 cm		> 120 cm	
cH (AT Maxi)	> 95 cm		> 105 cm		> 105 cm		> 125 cm		> 125 cm		> 135 cm	
cG siliconenband	43-48	49-57	45-52	53-62	49-56	57-67	53-60	61-72	56-64	65-77	60-68	69-81
cG	43 - 57		45 - 62		49 - 67		53 - 72		56 - 77		60 - 81	
cF	39 - 52		41 - 56		44 - 60		47 - 65		50 - 69		53 - 73	
cE	30 - 37		33 - 40		35 - 43		37 - 45		39 - 48		41 - 51	
cD	27 - 33		29 - 36		32 - 39		34 - 42		36 - 45		38 - 48	
cC	28 - 34		30 - 37		33 - 40		35 - 43		37 - 46		39 - 49	
cB1	23 - 27		24 - 29		26 - 32		29 - 35		31 - 37		33 - 39	
cB	18 - 20		20 - 22		22 - 24		24 - 26		26 - 28		28 - 30	
cY	26 - 31		28 - 33		29 - 35		31 - 37		32 - 38		33 - 40	
cA	17 - 22		19 - 24		21 - 26		23 - 29		25 - 32		27 - 34	

Longueur D

ℓD petite: 34 - 38 cm
ℓD regular: 38 - 45 cm

Longueur G

ℓG petite: 60 - 69 cm
ℓG regular: 69 - 80 cm
ℓG long*: 80 - 90 cm

Longueur K (panty)

ℓK petite: 65 - 75 cm
ℓK regular: 75 - 85 cm

Mesures en cm

*Seulement disponible en JOBST Classic